

Instructions for Use

Life Science Kits & Assays



smart DNA prep (m)

Order No.:

845-KS-8000010 10 reactions

845-KS-8000050 50 reactions

Publication No.: HB_KS-8000_e_1, \$&\$%

This documentation describes the state at the time of publishing.
It needs not necessarily agree with future versions. Subject to change!

Print-out and further use permitted with indication of source.

© Copyright 201, , Analytik Jena AG, AJ Innuscreen GmbH

Manufacturer:

AJ Innuscreen GmbH
Robert-Rössle-Straße 10
13125 Berlin
Made in Germany!

Distribution/Publisher:

Analytik Jena AG
Konrad-Zuse-Straße 1
07745 Jena · Germany

Phone +49 3641 77 9400
Fax +49 3641 77 767776
www.analytik-jena.com
info@analytik-jena.com

Contents

1	Introduction.....	3
1.1	Intended use.....	3
1.2	Notes on the use of this manual.....	4
2	Safety precautions.....	5
3	Storage conditions	6
4	Function testing and technical assistance.....	7
5	Product use and warranty.....	7
6	Kit components	8
6.1	Components not included in the kit.....	9
6.2	Components needed for isolation of nucleic acids from bacteria.....	9
6.3	Components needed for isolation of nucleic acids from yeasts	10
7	GHS Classification	11
7.1	Hazard phrases.....	12
7.2	Precaution phrases.....	13
8	Sample preparation for eukaryotic cells.....	14
8.1	Product specifications	14
8.2	Lysis of starting material.....	14
9	Sample preparation of tissue samples	15
9.1	Product specifications	15
9.2	Proteolytic lysis of starting material.....	15
10	Sample preparation of bacteria cell pellets.....	16
10.1	Product specifications	16
10.2	Resuspension of starting material.....	16
10.3	Pre-lysis of resuspended starting material.....	16
10.3.1	Gram-negative bacteria	16
10.3.2	Gram-positive bacteria.....	17

10.3.3 Staphylococcus	18
10.4 Proteolytic lysis step.....	19
11 Sample preparation of yeast cell pellets	20
11.1 Product specifications	20
11.2 Resuspension of starting material.....	20
11.3 Pre-lysis of resuspended starting material.....	20
11.4 Proteolytic lysis step.....	21
12 SmartExtraction protocol.....	22
12.1 Binding DNA to SE Macro Beads	22
12.2 Washing and removing of alcohol.....	22
12.3 Elution of DNA from SE Macro Beads	23
13 Troubleshooting.....	24
14 Related Products	25

1 Introduction

1.1 Intended use

The smart DNA prep (m) kit has been designed for the manual isolation of high molecular weight genomic DNA (200 kb - > 500 kb) from tissue samples, eukaryotic cells, rodent tails, bacteria and yeasts.

The procedure starts with the lysis of the starting material. After lysis the sample is transferred into a SmartExtraction Tube (SE Tube). The SE Tube contains a material with unique Smart Modified Surfaces which adsorbs the genomic DNA. After washing steps the nucleic acid is dissolved from the surface of the modified material and is now ready to use for downstream applications. The whole extraction process is simple to handle. The unique extraction chemistry in combination with Smart Modified Surfaces is optimized to get a maximum of yield and quality.









CONSULT INSTRUCTION FOR USE



This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of results cannot be guaranteed if there are any deviations from the instructions in this package insert.

1.2 Notes on the use of this manual

For easy reference and orientation, the manual uses the following warning and information symbols as well as the shown methodology:

Symbol	Information
	REF Catalogue number.
	Content Contains sufficient reagents for <N> tests.
	Storage conditions Store at room temperature, unless otherwise specified.
	Consult instructions for use This information must be observed to avoid improper use of the kit and the kit components.
	Expiry date
	Lot number The number of the kit charge.
	Manufactured by Contact information of manufacturer.
	For single use only Do not use components for a second time.
	Note / Attention Observe the notes marked in this way to avoid operating errors for obtaining correct results.

The following systematic approach is introduced in the manual:

- The chapters and figures are numbered consecutively.
- A cross reference is indicated with an arrow (e.g. → "Notes on the use of this manual" p. 4).
- Working steps are numbered.

2 Safety precautions

NOTE

Read through this chapter carefully prior to guarantee your own safety and a trouble-free operation.

Follow all the safety instructions explained in the manual, as well as all messages and information, which are shown.

All due care and attention should be exercised in handling the materials and reagents contained in the kit. Always wear gloves while handling these reagents and avoid any skin contact! In case of contact, flush eyes or skin with a large amount of water immediately.



FOR SINGLE USE ONLY!

This kit is made for single use only!

ATTENTION!

Don't eat or drink components of the kit!

The kit is designed to be handled only by educated personnel in a laboratory environment!

If the buffer bottles are damaged or leaking, wear gloves and protective goggles when discarding the bottles in order to avoid any injuries. This kit is to be used with potential infectious human samples. Therefore, all liquid waste must be considered as potentially infectious and must be handled and discarded according to local safety regulation.

Please observe the federal, state and local safety and environmental regulations. Follow the usual precautions for applications using extracted nucleic acids. All materials and reagents used for DNA should be free of DNases or RNases.

ATTENTION!

Do not add bleach or acidic components to the waste after sample preparation!

NOTE

Emergency medical information in English and German can be obtained 24 hours a day from:

Poison Information Center, Freiburg / Germany

Phone: +49 (0)761 19 240.

For more information, please ask for the material safety data sheets (MSDS's).

3 Storage conditions

Store lyophilized **Proteinase K** at 4 °C to 8 °C! Divide dissolved **Proteinase K** into aliquots and store at -22 °C to -18 °C. Repeated freezing and thawing will reduce the activity dramatically!

All other components of the smart DNA prep (m) kit should be stored dry at room temperature (15 °C to 30 °C). When stored at room temperature, the kit is stable until the expiration date printed on the label on the kit box.

If there are any precipitates within the provided solutions dissolve these precipitates by careful warming. Before every use make sure that all components have room temperature. For further information see chapter "Kit components" (→ p. 8).

4 Function testing and technical assistance

The Analytik Jena AG guarantees the correct function of the kit for applications as described in the manual. This product has been produced and tested in an ISO 13485 certified facility.

We reserve the right to change or modify our products to enhance their performance and design. If you have any questions or problems regarding any aspects of the smart DNA prep (m) kit or other Analytik Jena AG products, please do not hesitate to contact us. For technical support or further information in Germany please dial +49 36 41 / 77 94 00. For other countries please contact your local distributor.

5 Product use and warranty

The kit is not designed for the usage of other starting materials or other amounts of starting materials than those, referred to in the manual (→ "Intended use") (→ "Product specifications"). Since the performance characteristics of Analytik Jena AG kits have just been validated for the application described above, the user is responsible for the validation of the performance of Analytik Jena AG kits using other protocols than those described below. Analytik Jena AG kits may be used in clinical diagnostic laboratory systems after the laboratory has validated the complete diagnostic system as required by CLIA' 88 regulations in the U.S. or equivalent regulations required in other countries.

All products sold by the Analytik Jena AG are subjected to extensive quality control procedures and are warranted to perform as described when used correctly. Any problems should be reported immediately.

NOTE

For research use only!

6 Kit components

IMPORTANT

Store lyophilized **Proteinase K** at 4 °C to 8 °C! Divide dissolved **Proteinase K** into aliquots and store at -22 °C to -18 °C. Repeated freezing and thawing will reduce the activity dramatically!



STORAGE CONDITIONS

All other components are stored at room temperature.

	Σ 10	Σ 50
REF	845-KS-8000010	845-KS-8000050
SE Tube	10	50
Lysis Solution CBV	5 ml	30 ml
Proteinase K	For 2 x 0.3 ml working solution	For 2 x 1.5 ml working solution
Binding Optimizer	1 ml	3 x 1 ml
Washing Solution LS (conc.)	2 ml (final volume 10 ml)	12 ml (final volume 60 ml)
Elution Buffer	15 ml	2 x 30 ml
Manual	1	1
Initial steps	Washing Solution LS (conc.) Add 8 ml of 96-99.8 % ethanol to the bottle Washing Solution LS mix thoroughly and keep the bottle always firmly closed! Proteinase K Dissolve by addition of 0.3 ml of ddH ₂ O, mix thoroughly and store as described above.	Washing Solution LS (conc.) Add 48 ml of 96-99.8 % ethanol to the bottle Washing Solution LS mix thoroughly and keep the bottle always firmly closed! Proteinase K Dissolve by addition of 1.5 ml of ddH ₂ O, mix thoroughly and store as described above.

6.1 Components not included in the kit

- 1.5 ml and 2.0 ml tubes
- 96–98.8 % ethanol (molecular biology grade, undenaturated)
- 80 % ethanol (molecular biology grade, undenaturated)
- 2-Propanol (molecular biology grade)
- ddH₂O
- optional RNase A (10 mg/ml)
- magnetic rack (Analytik Jena AG, 845-MR-0600001)
- 1 x PBS Buffer (137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 1.8 mM KH₂PO₄)

6.2 Components needed for isolation of nucleic acids from bacteria

- Lysozyme (stock solution: 10 mg/ml (400 U/μl))
- Mutanolysin (stock solution: 0.4 U/μl)
- Lysostaphin (stock solution: 0.4 U/μl)
- TE-Buffer



Alternatively:

- innuPREP Bacteria Lysis Booster (Analytik Jena AG; 845-KA-1000050)

6.3 Components needed for isolation of nucleic acids from yeasts

- Yeast Digest Buffer
 - 50 mM potassium phosphate
 - 10 mM DTT
 - pH 7.5
- Lyticase (stock solution: 10 U/ μ l)

7 GHS Classification

Component	Hazard contents	GHS Symbol	Hazard phrases	Precaution phrases
Binding Optimizer	Acetic acid 10–≤25 %	 Warning	315, 319	101, 102, 103, 280, 305+351+338, 362, 302+352, 403+233, 501
Proteinase K	Proteinase, Tritirachium album serine	 Danger	315, 319, 334, 317, 335	101, 102, 103, 261, 280, 305+351+338, 342+311, 405, 501

7.1 Hazard phrases

- 315 Causes skin irritation.
- 317 May cause an allergic skin reaction.
- 319 Causes serious eye irritation.
- 334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
- 335 May cause respiratory irritation.

7.2 Precaution phrases

101	If medical advice is needed, have product container or label at hand.
102	Keep out of reach of children.
103	Read label before use.
234	Keep only in original container.
261	Avoid breathing dust/fume/gas/mist/vapors/spray.
280	Wear protective gloves/protective clothing/ eye protection/face protection.
362	Take off contaminated clothing.
405	Store locked up.
501	Dispose of contents/container in accordance with local/regional/national/international regulations.
302+352	IF ON SKIN: Wash with plenty of water.
342+311	If experiencing respiratory symptoms: Call a POISON CENTER/doctor.
403+233	Store in a well-ventilated place. Keep container tightly closed.
305+351+338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

8 Sample preparation for eukaryotic cells

8.1 Product specifications

Starting material:

- Eukaryotic cells (1×10^5 – 1×10^7)

8.2 Lysis of starting material

1. Collect the cells by centrifugation with parameters adequate for the cell type (e.g. 5 minutes at 2,500 x g) and discard the supernatant.
2. Add **150 µl PBS** to the cell pellet and resuspend the pellet as much as possible by intensive pipetting up and down.
3. Transfer the resuspended cells into the **SE Tube**.
4. Add **200 µl Lysis Solution CBV** and **40 µl Proteinase K**. Vortex shortly and incubate the **SE Tube** at 55 °C for 20 minutes in a thermal shaker continuously shaking with 1,200 rpm.

Lysis time of 20 minutes is often sufficient to get enough DNA, but the prolongation of lysis time to 1 hour is also possible.

NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

Proceed with "SmartExtraction protocol" on p. 22.

9 Sample preparation of tissue samples

9.1 Product specifications

Starting material:

- Tissue samples (1 mg–100 mg)
- Rodent tail (0.1 cm–1 cm)

9.2 Proteolytic lysis of starting material

1. Cut the starting material into small pieces and put it into a 1.5 ml reaction tube.
2. Add **400 µl Lysis Solution CBV** and **40 µl Proteinase K**. Vortex shortly and incubate the tube at 55 °C for 1 hour in a thermal shaker continuously shaking with 1,200 rpm.

Lysis time of 1 hour is often sufficient to get enough DNA, but the prolongation of lysis time to 3 hours is also possible.

3. After lysis centrifuge the tube at maximum speed for 5 minutes to pellet unlysed material. Carefully transfer the supernatant into the **SE Tube**.

NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

Proceed with "SmartExtraction protocol" on p. 22.

10 Sample preparation of bacteria cell pellets

10.1 Product specifications

Starting material:

- Bacteria cell pellets (1×10^5 – 1×10^9 cells)

10.2 Resuspension of starting material

1. Collect the cells by centrifugation with parameters adequate for the cell type (e.g. 10 minutes at 3,000 x g) and discard the supernatant.
2. Resuspend the bacteria cell pellet in **170 µl TE Buffer**. After resuspension start enzymatic pre-lysis as described below. Requirements for pre-lysis depend on the cell type.

10.3 Pre-lysis of resuspended starting material

10.3.1 Gram-negative bacteria

Although Gram-negative bacteria do not require a pre-lysis-step, using Lysozyme (not included in the kit) can enhance the efficiency of lysis.

Using Lysozyme

Stock solution of Lysozyme: 10 mg/ml (400 U/µl)

Add **20 µl Lysozyme** to the resuspended cells and incubate at 37 °C for 30 minutes under continuous shaking.

Proceed with “Proteolytic lysis step” on p. 19.

10.3.2 Gram-positive bacteria

Gram-positive bacteria require a pre-lysis-step using Mutanolysin and/or Lysozyme (not included in the kit).

Using Lysozyme

Stock solution of Lysozyme: 10 mg/ml (400 U/μl)

Add **20 μl Lysozyme** to the resuspended cells and incubate at 37 °C for 30 minutes under continuously shaking.

Using Mutanolysin

Stock solution of Mutanolysin: 0.4 U/μl

Add **5 μl Mutanolysin** to the resuspended cells and incubate at 37 °C for 30 minutes under continuously shaking.

Proceed with "Proteolytic lysis step" on p. 19.

NOTE

Lysozyme and Mutanolysin exert synergistic activity. Using both enzymes together will increase the yield of isolated nucleic acids.

Alternatively:

Use the **innuPREP Bacteria Lysis Booster** (Analytik Jena AG; 845-KA-1000050)

The innuPREP Bacteria Lysis Booster Kit has been developed for a highly efficient pre-lysis of bacterial cell walls by generating spheroblasts. This new mixture of different enzymes boost the lysis of all bacteria in particular the hard-to-lyse microorganisms like *Streptococcus*, *Lactobacillus*, *Staphylococcus*, *Bacillus* and *Clostridium*.

Prepare the enzyme mix according to the manual of the innuPREP Bacteria Lysis Booster.

Add **20 μl** of the prepared **enzyme mix** to the sample and vortex it shortly. Incubate the sample for 30 minutes at 37 °C.

Proceed with "Proteolytic lysis step" on p. 19.

10.3.3 Staphylococcus

For Lysis of *Staphylococcus* the enzyme Lysostaphin is recommended (not included in the kit)

Stock solution of Lysostaphin: 0.4 U/ μ l

Add **10 μ l Lysostaphin** to the resuspended cells and incubate at 37 °C for 30 minutes under continuously shaking.

Proceed with "Proteolytic lysis step" on p. 19.

Alternatively:

Use the **innuPREP Bacteria Lysis Booster** (Analytik Jena AG; 845-KA-1000050)

The innuPREP Bacteria Lysis Booster Kit has been developed for a high efficient pre-lysis of bacterial cell walls by generating sphaeroblasts. This new mixture of different enzymes boost the lysis of all bacterias in particular the hard-to-lyse microorganisms like *Streptococcus*, *Lactobacillus*, *Staphylococcus*, *Bacillus* and *Clostridium*.

Prepare the enzyme mix according to the manual of the innuPREP Bacteria Lysis Booster.

Add **20 μ l** of the prepared enzyme mix to the sample and vortex it shortly. Incubate the sample for 30 minutes at 37°C.

Proceed with "Proteolytic lysis step" on p. 19.

10.4 Proteolytic lysis step

Transfer the pre-lysed cells into the **SE Tube** and add **200 µl Lysis Solution CBV** and **30 µl Proteinase K**. Vortex shortly and incubate the SE Tube at 55 °C for 30 minutes in the a thermal shaker continuously shaking with 1,200 rpm.

Lysis time of 30 minutes is often sufficient to get enough DNA. If the sample is not clear after 30 minutes prolong the incubation time until the sample is clear.

NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

11 Sample preparation of yeast cell pellets

11.1 Product specifications

Starting material:

- Yeast cell pellets (1×10^5 – 1×10^9 cells)

11.2 Resuspension of starting material

1. Collect the cells by centrifugation with parameters adequate for the cell type (e.g. 10 minutes with 3,000 x g) and discard the supernatant.
2. Resuspend the yeast cell pellet in **200 µl Yeast Digest Buffer** (→"Components needed for isolation of nucleic acids from yeasts" p. 10). After resuspension start enzymatic pre-lysis as described below.

11.3 Pre-lysis of resuspended starting material

For lysis of yeast cells the enzyme Lyticase is recommended (not included in the kit)

Stock solution of Lyticase: 10 U/µl

Add **10 µl Lyticase** (not included in the kit) to the resuspended cells and incubate at 37 °C for 30 minutes under continuous shaking.

Proceed with "Proteolytic lysis step" on p. 21.

11.4 Proteolytic lysis step

Transfer the pre-lysed cells into the **SE Tube** and add **200 µl Lysis Solution CBV** and **30 µl Proteinase K**. Vortex shortly and incubate the **SE Tube** at 55 °C for 30 minutes in a thermal shaker continuously shaking with 1,200 rpm.

Lysis time of 30 minutes often is sufficient to get enough DNA. If the sample is not clear after 30 minutes prolong the incubation time until the sample is clear.

NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

12 SmartExtraction protocol

12.1 Binding DNA to SE Macro Beads

1. After lysis add **40 µl Binding Optimizer** and **350 µl 2-Propanol** to the sample in the SE Tube.
2. Place the **SE Tube** into a thermal shaker and incubate for 3 minutes at 1,400–1,800 rpm.
3. Place the **SE Tube** into a magnetic rack for separation of the SE Macro Beads. Discard the supernatant.

12.2 Washing and removing of alcohol

1. Add **800 µl of Washing Solution LS** and wash the SE Macro Beads by inverting the **SE Tube within the magnetic rack** five times. Discard the Washing Solution LS (do not remove the SE Tube from the magnetic rack).
2. Add **800 µl of 80 % ethanol** and wash the SE Macro Beads by inverting the **SE Tube with the magnetic rack** five times. Discard the ethanol (do not remove the SE Tube from the magnetic rack).
3. Add **800 µl of 80 % ethanol** and wash the SE Macro Beads by inverting the **SE Tube within the magnetic rack** five times. Discard the ethanol (do not remove the SE Tube from the magnetic rack).

4. Let stay the **SE Tube** in the magnetic rack and wait some seconds. Remove residual ethanol as much as possible using a pipet (also remove ethanol which can be inside the lid).
5. Place the **SE Tube** with opened lid in a thermal shaker and incubate for 15 minutes at 65 °C shaking with 400 rpm to remove the ethanol completely. If the ethanol is not completely removed prolong the incubation time.

12.3 Elution of DNA from SE Macro Beads

1. Add 200 µl–1,000 µl of **Elution Buffer** to the **SE Tube**. The amount depends on starting amount of sample and expected yield. Incubate the **SE Tube** in a thermal shaker for 15 minutes at 65 °C shaking with 1,000 rpm. If the DNA is not completely dissolved (DNA is visible as a clump or a strand) a prolongation of the elution step is strongly recommended.
2. Place the **SE Tube** into the magnetic rack and transfer the DNA into a new tube.

IMPORTANT NOTE

The extracted DNA is of high molecular weight (200 kb–500 kb). Therefore, the DNA might be very viscous. The dissolving step is crucial for successful extraction and for a maximum of yield. If the DNA content is too high, increase the amount of Elution Buffer and prolong the elution step.

If you do not need high molecular weight DNA you can shear the DNA e.g. by using ultrasound or by passing the eluate through a needle.

13 Troubleshooting

Problem / probable cause	Comments and suggestions
Low amount of extracted DNA	
Insufficient lysis	Increase lysis time. Reduce amount of starting material.
Incomplete elution	Prolong the incubation time with Elution Buffer.
Preparation without Binding Optimizer	It is important to add the Binding Optimizer to the lysed sample as described. Binding Optimizer need to be added after lysis of sample is finished!
High viscosity extracted DNA	
Insufficient amount of Elution Buffer	Elute the DNA with a higher volume of Elution Buffer.
Degraded or sheared DNA	
Old material	Old material often contains degraded DNA.
Pipetting steps performed too rigorous	Pipet more carefully and/or use wide bore pipette tips.

14 Related Products

Name	Amount	Order No.
Additional products for nucleic acid purification		
innuPREP Bacteria Lysis Booster	50 rxn	845-KA-1000050
innuPREP Proteinase K	6 mg	845-CH-0010006
	30 mg	845-CH-0010030
Kits for Nucleic acid purification		
smart Blood DNA Midi prep (m)	10 rxn	845-KS-8100010
	50 rxn	845-KS-8100050
Equipment		
Magnetrack 6-fach	For up to six samples	845-MR-0600001
TS1 ThermoShaker, 100–240 V	EU plug	846-051-500
	US plug	846-051-590
Block Module 24 x 2,0 ml tubes for TS1		846-051-516

Headquarters

Analytik Jena AG
Konrad-Zuse-Str. 1
07745 Jena · Germany

Phone +49 3641 77 70
Fax +49 3641 77 9279
info@analytik-jena.com
www.analytik-jena.com

Pictures: Analytik Jena AG
Subject to changes in design and scope of delivery as well as further technical development.